


For a better hearing experience



The Pockettio helps you to correct hearing impairments and to better understand your environment. Whether you are with family or friends, listening to your MP3 player or watching television – the Siemens Pockettio will improve your hearing in every situation. This easy to use hearing aid will allow you to, once again, actively take part in day-to-day life.

Local Contact Information

Siemens contact
House number and street name
City and postal code
Country
Phone: +XX-XXX-XXX-XXXX
www.XXXXXXXXXXXXXXXXXX

Legal Manufacturer

Siemens Audiologische Technik GmbH
Gebbertstrasse 125
DE-91058 Erlangen
Germany

Global Business Unit

Siemens Audiologische Technik GmbH
Gebbertstrasse 125
91058 Erlangen
Germany
Phone: +49 9131 308 0

Global Siemens Headquarters

Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens Healthcare Headquarters

Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany
Phone: +49 9131 84 0
www.siemens.com/healthcare

www.siemens.com/hearing



Pockettio Digital

Rediscover the feeling of being able to hear – with the Siemens Pockettio.

Answers for life.

SIEMENS

Pockettio Digital


Equipped for all events

The Pockettio from Siemens is a robust, digital hearing system that is simple and convenient to use. The user can choose between three programs depending on the hearing situation. The device features modern technology that a hearing care professional can adjust to accommodate each individual's hearing loss. The Pockettio boasts top-notch Siemens quality at an outstanding price.

The Siemens brand has been a symbol of German inventiveness for more than 130 years. Siemens Audiology is committed to developing solutions for better hearing. In order to ensure best possible hearing for everyone, Siemens works together with universities, clinics and hearing centers.



All details at a glance



- | | |
|--|---|
| 1 ON/OFF switch | 5 Audio input, socket to connect to TV or other audio devices |
| 2 Volume control with clear marking on top and front | 6 Program change button |
| 3 Battery door | 7 Receiver cord |
| 4 LED indicator | 8 Receiver with ear tip |

Optional accessories (not included in initial packaging):

- Audio input cable, which connects the Siemens Pockettio to TV, MP3 player and other audio devices
- Y-cable for binaural listening

Easy to wear

- A 3 receiver tips for comfortable fit – suits every ear canal
- B Lanyard for convenient carrying
- C Clip for fastening on belt or pocket



Easy to use

- D Individual change of loudness with volume control
- E Optical (LED indicator) and acoustical signals for power on, program selection, low battery capacity
- F Easy to operate battery door and battery insertion, long battery life and use of standard battery (AA)



- G 3 programs for different hearing situations:
 - Program 1 optimized for speech understanding
 - Program 2 optimized for noisy surroundings
 - Program 3 optimized for hearing comfort



- H Audio input for connection to TV or other audio devices, automatically switches to TV audio when Siemens Pockettio is connected to TV via cable