



Be one with your world.

Introducing Siemens Tek. The wireless enhancement for Siemens hearing instruments.

www.siemens.com/hearing

SIEMENS



Your music.

Your movies.

Your family.

Your friends.

Your shows.

Be one.

Your info.

Your world.



You live in a connected world. Mobile. Instant. Always on. Every day there's some new way to be entertained, be informed, be in touch. And now, with Siemens Tek™, there's a way to bring it wirelessly together, with your hearing instruments, for the richest listening experience possible.



TEK)))

Look for this symbol
for Tek-ready products.


Now it's so easy to
link to your
world.

Tek Connect is your key to controlling wireless communications between your hearing instruments and your cell phone, your audio system, your TV, even your PC.

Tek Connect is the remote control for Siemens Tek, the wireless enhancement for select Siemens hearing instruments. Specifically designed to meet the demands of today's interactive hearing instrument wearer, Tek Connect speaks the same language as all of your other high-tech toys. Giving you better performance and sound quality. Not to mention convenience and comfort. Make and take hands-free cell phone calls directly through your hearing instruments. Listen to your MP3 files in natural, stereophonic, digital sound. Watch TV and DVDs with optimized audio. Tek Connect. Wireless. Effortless. Seamless. It's time for everything to get in sync.

Link to your mobile.

Now there's an easier way to connect with your mobile phone. Answer calls from a button on Tek Connect. Get the signal directly in both hearing instruments. Leave your phone in your pocket. Life is easy.




1. Tek Connect communicates wirelessly with your hearing instruments, turning them into a virtual hands-free headset.
2. Tek Connect tells you when a call comes in, even when you're in a different program mode.
3. Make, take, and end calls directly from your Tek Connect to your Bluetooth-enabled mobile or cordless phone.
4. Link to one phone at a time, and enjoy clear, digital, disturbance-free calls. Simple. Easy. Convenient.





Link to your television.


Enjoy optimum audio whether you're screening a DVD or simply watching TV. There's no disturbance. No feedback. With sound streamed directly to your hearing instruments, there's no need to fight over the volume control. Sounds good, doesn't it?



1. For all practical purposes, Tek Connect turns your hearing instruments into perfectly balanced high-tech, stereo headphones.
2. Tek Connect wirelessly transmits stereo audio signals simultaneously to both hearing instruments.
3. Tek Connect communicates with the Tek Transmitter using a proprietary Bluetooth protocol developed by Siemens. With ordinary Bluetooth, there is an inherent transmission delay that produces an annoying echo, and can cause the picture and sound to be slightly off. But with Siemens Tek, signal delay is reduced so your eyes and ears are always in sync.
4. The Tek Transmitter is connected to your TV and acts as a high-performance base station with a range of 10 meters (about 40 feet). It comes pre-paired with Tek Connect, so simply turn it on. Sit back. Relax. Enjoy the show.

Link to your music.

From now on, nothing will get between you and your music. Tek Connect streams stereo sound wirelessly from your audio system or MP3 player right to your hearing instruments. Suddenly earbuds are so 20th-century.



1. Tek Connect transforms your hearing instruments into a precision-paired, high-performance, stereo headset.
2. Using a proprietary transmission protocol, Tek Connect streams your music to both hearing instruments in the same instant.
3. Your MP3 player or stereo can be connected directly to Tek Connect using the 2.5 mm stereo input plug.
4. Tek Connect can also communicate with your MP3 player or stereo through the Tek Transmitter. It has a signal radius of 10 meters (about 40 feet). In case you feel like dancing.





Your life in

stereo.



If the sounds of your life are in stereo, shouldn't you hear them that way? Now you can. With Siemens Tek.

The technology of Siemens Tek helps your hearing instruments deliver true, natural hearing performance in any listening environment. So what you hear sounds more natural. In stereo. In your hearing instruments. All in a microsecond. Instantly.

Only Siemens Tek gives you control of the high and low frequencies of your hearing instruments. Adjust the timbre of the sound you hear directly from your Tek Connect.

Fine-tune the treble and bass ranges with the touch of a fingertip. And enrich the overall quality of your personal sound.

Siemens Tek automatically synchronizes the high-definition, balanced, binaural signals between your instruments. No need for cables or headphones. No interference from other sounds or noises around you. The result is a hearing experience controlled by you, for you. Plug. Play. Perfect.



Isn't it nice
when everything
plays
together?

These days, to make the best technology you have to understand you're not the only technology. Siemens Tek links you to all your favorite devices, so you can be connected to all of your favorite things. With a 130-year history in audiology, Siemens is committed to giving you the kinds of hearing enhancements that match the way you live. So you can truly be one with your world.

Research In Motion, the RIM logo, BlackBerry, the BlackBerry logo and SureType are registered with the U.S. Patent and Trademark Office and may be pending or registered in other countries – these and other marks of Research In Motion are used under license.

The information in this document contains general descriptions of the technical options available, which do not always have to be present in individual cases and are subject to change without prior notice.

The required features should therefore be specified in each individual case at the time of conclusion of the respective contract.

© 02.2008 Siemens AG
Order No. A91SAT-00396-99C1-7600
Printed in Germany
104300/10054 02085.

Siemens Audiologische Technik GmbH
Gebbertstrasse 125
91058 Erlangen
Germany
Phone +49 9131 308-0

www.siemens.com/hearing